

DISPOSABLE MASK TYPE IIR



Disposable facial protection mask with 3 layers of protection and soft fixation on the ears by means of elastics, with integrated nasal support inside (easy-adjust nose clip). Blue exterior and white interior.

Non-sterile product. Single use only. Not reusable.

Boxes of 50 units

- ✓ 3 layers
- ✓ Type IIR
- ✓ With adjustable nasal clip
- ✓ Soft elastics



Declaration of Conformity available at www.chemitool.com

TECHNICAL SPECIFICATIONS

Model	Disposable mask type IIR
Sizes	One size
Quantity per dispenser	50 units
Packaging	40 dispensers 2 000 máscaras
Material	3 layers: 20gr blue TNT, Meltblown and 20gr white TNT.
Mask Size	17.5 cm x 9.5 cm

INSTRUCTIONS FOR USE



1. Hold the mask with the blue side facing out.
2. Flip the clip area up and place it against the bridge of your nose.
3. Attach the elastics to the ears.
4. Make sure the mask completely covers your nose and mouth.
5. Press the nose clip to fit your nose.

BEFORE USING THE MASK, WASH YOUR HANDS AND MAKE SURE THEY ARE NOT DAMAGED.

STORAGE

Store out of direct sunlight, away from chemicals, and in a dry environment.

SAP	TAMANHO	EAN
CHS07050404	one size	5604630061773

Note: Technical information provided, whether orally or in writing, is based on our current knowledge and should be considered collaboration without obligation. The use of the product is beyond our control, therefore, we exclude all responsibility for the improper use of the product. The customer is responsible for confirming and evaluating (through testing) whether the product is suitable for the process and type of use in question. Our purpose is exclusively to guarantee the quality of the products, in accordance with our standards.